



2. Fix in 1ml cold 70% ethanol. Add drop wise to cell pellet while vortexing. This should ensure fixation of all cells and minimise clumping.
3. Fix for at least 30 minutes on ice. Specimens can be left at this stage for several weeks (make sure you seal the tubes for long term storage).
4. Pellet cells at higher speed for 5 minutes; decant the supernatant being careful not to lose the pellet. Note that ethanol-fixed cells require higher centrifugal speeds to pellet compared to unfixed cells since they become more buoyant upon fixation.
5. Wash once with PBS + 0.1% Triton X solution.
6. Wash once with PBS + 0.1% Triton X + 3% BSA.
7. Pour off supernatant and add directly to the pellet 2 $\mu$ l of rabbit Anti-Phospho Histone 3 antibody, mix well and incubate at room temperature for 45 minutes.
8. Wash twice with PBS + 0.1% Triton X solution.